

Select the right dressing based on wound conditions

Product Name

Exudate Level

Indications

Protection & Fixation



Klario Film
Transparent Film Dressing
(Frame Delivery)



Wound protection,
fixation

Post-Op Care



Klario Post-Op
Transparent Waterproof
Dressing with Absorbent
Pad



Post-operative
incisions, cuts, minor
abrasions



Regeneo Hydro Fix
Hydrocolloid Dressing
with Transparent Film
Border



Superficial wounds,
pressure ulcers, burns,
post-operative
wounds

Moisture & Exudate Management



Regeneo Hydro
Hydrocolloid Dressing



Pressure ulcers,
burns, superficial
wounds

Pilih balutan yang tepat sesuai dengan kondisi luka

Product Name	Exudate Level	Indications
Moisture & Exudate Management		
	Regeneo Foam Non-Adhesive Foam Dressing Pad	 Moderate exudate wounds, pressure ulcers, venous ulcers
	Regeneo Foam Non-Adhesive Foam Dressing Pad with Film Layer	 Moderate exudate wounds, pressure & venous ulcers, post-operative wounds
	Regeneo Alginate Calcium Alginate Dressing	 Minor bleeding, deep wound packing / cavity filling